Gynecare

TVT ABBREVO™

Continence System

GYNECARE TVT ABBREVO™ kontinenssystem

**GYNECARE TVT ABBREVO™-
continentiesysteem**

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

CE 2797

P24070
LAB0012890v8
01/2021
Made in Switzerland
© Ethicon, Inc. 2021


Ethicon, SÀRL
Puits Godet 20
Neuchâtel
CH-2000
Switzerland
1-877-ETHICON
+1-513-337-6928


1


2


3


4


5


6


7


8A


8B


9


10


11


12A

12B

13

14

Continence System

GYNECARE TVT ABBREVO™ Continence System

GYNECARE TVT ABBREVO™ Implant Assembly
Sterile, Single-Patient Use

GYNECARE TVT ABBREVO™ Placement Loop
Sterile, Single-Patient Use

GYNECARE TVT ABBREVO™ Helical Passers
Sterile, Single-Patient Use

GYNECARE TVT ABBREVO™ Atraumatic Winged Guide
Sterile, Single-Patient Use

Please read all information carefully.

Failure to properly follow instructions may result in improper functioning of the device and/or may lead to injury. Properly following instructions will not eliminate the risk of adverse events.

Introduction

This package insert is designed to provide instructions for use of the GYNECARE TVT ABBREVO Continence System, including the GYNECARE TVT ABBREVO Device, Helical Passers and Atraumatic Winged Guide. It is not a comprehensive reference to surgical technique for correcting Stress Urinary Incontinence (SUI). This is not a technique manual nor a substitute for appropriate training and experience in surgical technique for correcting SUI. The device is intended to be used only by surgeons trained in the surgical treatment of SUI, and specifically in the use of the GYNECARE TVT ABBREVO Continence System. Training should be in line with the training guidelines specified by the International Urogynecological Association (IUGA) or by other international, national, or local Institutions. Adverse reactions may occur notwithstanding the training and experience of the surgeon. These instructions are intended for general use of the GYNECARE TVT ABBREVO Continence System. The removal of the product from the patients (if needed) should be performed by a trained urogynecologist, urologist or gynecologist experienced in the removal of pelvic floor meshes. Variations in use may occur in specific procedures due to individual technique and patient anatomy.

Description

The GYNECARE TVT ABBREVO Continence System is a sterile, single-patient use procedure kit consisting of:

GYNECARE TVT ABBREVO Implant Assembly

The GYNECARE TVT ABBREVO Implant Assembly consists of one piece of blue (up to 0.28 weight percent phthalocyanine blue, Color Index number 74160) PROLENE™ Polypropylene Mesh (Implant), approximately ½ inch x 4.7 inches (1.1 cm x 12 cm), covered by clear plastic Sheaths and held between two Helical Passer Sheaths (white polyethylene tube receptacles). Based upon the maximum implanted mesh length, the maximum amount of dye (phthalocyanine blue, Color Index number 74160) is 0.00037 g. The Helical Passer Sheaths are bonded to the mesh Implant sheaths and the Positioning Lines (PROLENE™ Polypropylene Monofilament). PROLENE Mesh is constructed of knitted monofilaments of extruded polypropylene, identical in composition to that used in PROLENE™ Polypropylene Nonabsorbable Surgical Sutures. The Implant is approximately 0.027 inches (0.7 mm) thick. This material, when used as a suture, has been reported to be non-reactive and to retain its strength indefinitely in clinical use. PROLENE Mesh is knitted by a process which interlinks each fiber junction.

GYNECARE TVT ABBREVO Placement Loop

The GYNECARE TVT ABBREVO Placement Loop is a sterile, single-patient use device consisting of a loop of PROLENE Polypropylene Monofilament with an attached polypropylene button. This loop and button are pre-assembled as part of the GYNECARE TVT ABBREVO Implant Assembly at the center of the mesh to aid in symmetric placement of the mesh.

GYNECARE TVT ABBREVO Helical Passers

The GYNECARE TVT ABBREVO Helical Passers are stainless steel, curved wire passers with plastic handles that are designed to deliver the GYNECARE TVT ABBREVO Implant. The Helical Passers are provided as left and right units, pre-assembled to the GYNECARE TVT ABBREVO Implant Assembly. The Helical Passer MUST NOT be bent or deformed in any way.

GYNECARE TVT ABBREVO Atraumatic Winged Guide

The GYNECARE TVT ABBREVO Atraumatic Winged Guide is a stainless steel accessory instrument which facilitates consistent passage of the GYNECARE TVT ABBREVO Helical Passers through the dissection tract. The Winged Guide is marked with an Insertion Zone to aid the surgeon's assessment of inserted depth. The Insertion Zone indicates a distance of 3 cm to 4 cm from the tip of the Winged Guide.

Indications / Intended Use

The GYNECARE TVT ABBREVO Continence System is intended for use in women as a sub-urethral sling for the treatment of SUI resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

The GYNECARE TVT ABBREVO Helical Passers and Atraumatic Winged Guide are intended to aid in the placement of the GYNECARE TVT ABBREVO Device.

Contraindications

As with any suspension surgery, this procedure should not be performed in pregnant patients. Additionally, because the PROLENE Polypropylene Mesh will not stretch significantly, this procedure should not be performed in patients with future growth potential, including women with plans for future pregnancy.

Patient Target Group and Factors

GYNECARE TVT ABBREVO is intended for use for non-pregnant adult female patients not planning future pregnancies, who are affected by SUI resulting from urethral hypermobility and/or intrinsic sphincter deficiency. Surgeons should use their surgical experience and judgment to determine if PROLENE Mesh is appropriate for certain patients. Patient-specific factors may impair wound healing, which may increase the likelihood of adverse reactions.

Application / Instructions for Use

Note: *Hand positions used in surgery may vary from those shown in the illustrations.*

1. Place the patient in the lithotomy position with the hips hyperflexed over the abdomen. The tip of the coccyx should be positioned flush with the edge of the table. Abduction of the legs should be performed after flexion of the hips is realized (see Figure 1).
2. The procedure can be performed under local, regional, or general anesthesia.
3. If desired, retract the labia to provide additional exposure.
4. Insert a urethral catheter into the bladder and empty the bladder.
5. If desired, mark the exit points of the device by tracing a horizontal line at the level of the urethral meatus, and a second line parallel to and 2 cm above the first line. Locate and mark the exit points on this second line 2 cm lateral to the folds of the thigh (the skin may be flattened by stretching) (see Figure 2).
6. Using Allis clamps for traction, make a 1 cm midline incision in the vaginal mucosa, starting 1 cm proximal to the urethral meatus, and two lateral sharp dissections.

Note: *It is suggested that device insertion be completed on one side before beginning the scissor dissection of the second side.*

7. After initiating sharp dissection, continue making a narrow blunt dissection toward, but NOT into, the obturator membrane using a “push-spread technique” with fine-pointed, curved scissors. The path of the lateral dissection should be in the horizontal plane, directed toward the ischio-pubic ramus in a 45° angle in relation to the coronal plane. Continue dissection just past the ischio-pubic ramus and into the obturator internus muscle, but do not perforate the obturator membrane. The aim of the dissection is to create a channel that is just wide enough for insertion of the Winged Guide (see Figure 3).
 - 3A. Horizontal Plane
 - 3B. Ischio-Pubic Ramus
 - 3C. Midline Plane
 - 3D. Vaginal Dissection

Important Note: *Dissection of the obturator membrane or a dissection channel that is wider than described may result in a reduction in fixation force.*

8. Remove the internal package workstation from the external package. Then remove the Winged Guide from the package workstation (see Figure 4).
 - 4A. Winged Guide
 - 4B. Placement Loop
 - 4C. Helical Passer/Implant Assembly
9. Insert the Winged Guide into the dissected tract next to the ischio-pubic ramus and into the obturator internus muscle, **but do not perforate the obturator membrane**. The Insertion Zone on the Winged Guide should be visible during this part of the procedure. Placement beyond the Insertion Zone may allow unintended entry into the obturator membrane or space of Retzius. If bony contact is not achieved after insertion of the Guide within the Insertion Zone, remove and re-evaluate the angle of dissection. If difficulty is encountered during insertion of the Guide, reconfirm the direction of the dissected tract with the scissors.

Note: *The open side of the Winged Guide must be facing the surgeon and be pressed tightly against the patient's contra-lateral thigh (see Figure 5).*

10. The Winged Guide should be held in place by the surgeon or assistant.
11. Remove the Helical Passer/Implant Assembly from the workstation (see Figure 4).
 - 4A. Winged Guide
 - 4B. Placement Loop
 - 4C. Helical Passer/Implant Assembly

Note: *To ensure correct orientation of the system, verify that the GYNECARE logo and thumb indent on the Helical Passer handles are facing the surgeon, and that the tips of the Helical Passer Sheaths are pointing outward. The Helical Passer in the surgeon's left hand must be used on the patient's right side. (see Figure 6)*

12. Place one of the Helical Passers on the sterile drape or other suitable sterile location until needed.
13. Insert the first Helical Passer into the dissected tract following the channel of the Winged Guide. Ensure that the Helical Passer handle is oriented so that the straight tip of the Helical Passer is aligned with the channel in the Guide and remains in that orientation until the tip reaches/contacts

the superior border of the inferior ischio-pubic ramus and is pushed slightly into the internus muscle, such that the tip of the Helical Passer is not into the obturator membrane.

14. Once Helical Passer/Implant Assembly is in position (i.e., when the tip of the Helical Passer is firmly positioned behind the ischio-pubic ramus), maintain handle position while advancing the tip and radius of the Helical Passer to perforate the obturator membrane (determined by a palpable pop), then remove the Winged Guide and retain it in the sterile field for later use in the procedure.
15. Once the Winged Guide has been removed, continue to slowly rotate the Helical Passer in an antero-medial direction as the Helical Passer is brought near the patient's midline, hugging the ischio-pubic ramus. Rotation in this manner while ensuring continuous contact with the bone allows a tight passage around the ramus (see Figures 7 and 8).
 - 7A. Starting position
 - 7B. Ending position
 - 8A. Starting position, side view
 - 8B. Ending position, side view

Note: *Never allow the handle to be oriented horizontal to the floor.*

16. The device should exit near or infero-medial to the previously determined exit points (if marked, see Figure 9). Making a skin incision at the site where the Helical Passer tents the skin may facilitate exiting of the device. When the tip of the Helical Passer Sheath appears at the skin opening, grasp its pointed tip with a clamp and, while stabilizing the Helical Passer Sheath near the urethra with the thumb, remove the Helical Passer by initially pushing the Helical Passer Sheath cranially with the thumb, then reversing the rotation of the handle. (see Figure 10)
17. Pull the Helical Passer Sheath completely through the skin until the Positioning Line and mesh sheath appear, leaving between 5 cm and 10 cm of Implant Assembly outside the vaginal incision before the Placement Loop. (see Figure 11)
18. Repeat the technique on the patient's other side, ensuring that the mesh will lie flat under the urethra and the Placement Loop is central between equal amounts of exposed mesh and mesh sheath. If cystoscopy was performed following the passage of the first Helical Passer, make sure the bladder is emptied prior to initiating passage on the second side. (see Figure 12A)

Note: *At this point cystoscopy may be performed at the discretion of the surgeon.*

Note: *If a twist in the mesh implant is discovered, ensure that the twist and Placement Loop are not positioned against the urethra as the mesh sheath is being pulled into position.*

19. When both Helical Passer Sheaths have been extracted through the skin incisions, the implant should be further positioned as follows: (see Figure 12)
 - a. Advance both mesh sheaths until the Placement Loop is central to the vaginal incision and the mesh is in contact with the urethra. (see Figure 12B)

Note: *Premature removal of the mesh sheaths may make subsequent adjustments difficult.*

- b. Cut the mesh sheath and Positioning Lines just medial to the Helical Passer Sheath on each side of the device beyond the skin exit points.

Note: *It is suggested to use a suitable means of counter-traction during sheath removal to ensure that the Placement Loop stays central. (This may reduce the risk of unintended over-tensioning of the mesh implant during sheath removal.)*

Note: *The Placement Loop should not be used as a means of counter-traction. Pulling on the Placement Loop to adjust mesh tension may damage the mesh and should not be done.*

- c. While maintaining suitable counter-traction of the mesh under the urethra, remove both sheaths by pulling them laterally. Take care that you DO NOT remove the two sets of Positioning Lines at this time. After sheath removal, the mesh should be in contact with the urethra once the suitable means of counter-traction are removed. Appropriate final positioning is at the discretion of the surgeon.
 - d. If the mesh implant needs to be loosened, the use of a blunt instrument is advisable, as the Placement Loop is not intended to be used to adjust the tensioning of the mesh.

Note: *Pulling on the Placement Loop to adjust mesh tension may damage the mesh and should not be done.*

- e. If the mesh requires additional tensioning, determine the best side to adjust and pull on BOTH ends of the Positioning Lines on that side to position the mesh appropriately while maintaining mesh symmetry during adjustment. (see Figure 13)

Note: *Failure to pull BOTH ends of the lines will detach the Positioning Lines from the mesh implant and no longer allow for mesh adjustment on that side.*

Note: *At this stage, it is advisable to perform a cough test or Credé maneuver (i.e., exerting manual pressure on the suprapubic area, in particular if the patient is under general anesthesia) after instilling the bladder with 300 mL (or up to the patient's preoperatively determined bladder capacity) to adjust mesh implant tensioning to the individual patient's characteristics.*

- f. When satisfied with mesh placement, remove the Positioning Lines from both skin exit points by gently pulling on ONE end of each Positioning Line.
 - g. The Placement Loop should now be removed by cutting the loop and removing both the loop and the button from the mesh. (see Figure 14)
20. Following mesh implant adjustment and removal of the Placement Loop, close the vaginal incision. Close the skin incisions with suture or surgical skin adhesive.

21. Post-operative indwelling catheterization is not typically required, but the bladder should be emptied at the end of the procedure. The patient should be encouraged to try to empty the bladder 2 to 3 hours after the procedure.

Disposal

Discard unintentionally opened/ partially used/ used devices and packaging according to your facility's policies and standard operating procedures concerning biohazardous materials and waste.

Warnings and Precautions

- Do not use the GYNECARE TVT ABBREVO Continence System in patients who are on anti-coagulation therapy.
- Do not use the GYNECARE TVT ABBREVO Continence System in a patient who has a urinary tract infection.
- Surgeons should be familiar with surgical technique for mid-urethral sling placement and should be adequately trained in the GYNECARE TVT ABBREVO Device procedure before employing the GYNECARE TVT ABBREVO Device. It is important that the Implant be located without tension under mid-urethra.
- Acceptable surgical practice should be followed for the procedure as well as for the management of contaminated or infected wounds.
- The procedure should be performed with care to avoid large vessels, nerves, bladder, and bowel. Attention to patient anatomy and correct passage of the device may minimize certain risks.
- Bleeding may occur post-operatively. Observe for any symptoms or signs before releasing the patient from hospital.
- Although bladder injury is unlikely to occur with this technique, cystoscopy may be performed at the discretion of the surgeon.
- Do not remove the plastic sheaths until the Implant has been properly positioned.
- Ensure that the mesh implant is positioned symmetrically and in contact with the mid-urethra with minimal tension to achieve desired results.
- Do not perform this procedure if you think the surgical site may be infected or contaminated.
- Since no clinical information is available about pregnancy following sub-urethral sling procedure with the GYNECARE TVT ABBREVO Device, the patient should be counselled that future pregnancies may negate the effects of the surgical procedure and the patient may again become incontinent.
- Since limited clinical information is available with delivery following a sub-urethral sling procedure with the GYNECARE TVT ABBREVO Device, in case of pregnancy, mode of delivery should be determined by the obstetrician in consultation with the surgeon.
- Post-operatively, the patient should be advised to rest for the first 24 to 48 hours. Further, the patient should be advised to refrain from heavy lifting and/or exercise (e.g., cycling, jogging) for at least four to six weeks and from intercourse for at least one month. The patient can usually return to other normal activity after two weeks.
- Should dysuria, bleeding, or other problems occur, the patient should be instructed to contact the surgeon immediately.
- Transient leg pain lasting 24 to 48 hours may occasionally occur and can usually be managed with mild analgesics.
- As with other incontinence procedures, de novo detrusor instability may occur following a sub-urethral sling procedure utilizing the GYNECARE TVT ABBREVO Device. To minimize this risk, make sure to place the mesh implant as described above.
- Do not contact the PROLENE Mesh with any staples, clips or clamps as mechanical damage to the mesh may occur.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and surgeons.
- Prophylactic antibiotics can be administered according to the surgeon's usual practice.
- In the event of a damaged product, the product should be discarded and a new one retrieved to complete the procedure. In the event of a product malfunction during use, continue or discontinue use at surgeon's discretion.

Adverse Reactions / Undesirable Side Effects

- Punctures or lacerations of vessels, nerves, structures, or organs, including the bladder, urethra, or bowel, may occur and may require surgical repair.
- Local irritation at the wound site may occur.
- As with any implant, a foreign body response will occur, the extent of which may differ. This response could result in extrusion, erosion, exposure, fistula formation and/or chronic inflammation, the severity of which is unpredictable, or other adverse reactions, which may be ongoing.
- Fistula formation, acute and chronic inflammation and ongoing risk of mesh extrusion, exposure, or erosion into the vagina or other structures or organs (such as bladder, urethra or rectum) which may be difficult to treat and result in consequent pain. Mesh extrusion, exposure, or erosion into the vagina may also cause offensive vaginal discharge.
- Infection following transvaginal implantation. As with all surgical procedures and the implantation of foreign bodies, there is a risk of infection and PROLENE Mesh may potentiate an existing infection.
- Pain - which may be severe and chronic
- Temporary or chronic voiding dysfunction (or difficulty voiding) or urinary retention/obstruction independent from that caused by overcorrection or urethral hypermobility, i.e., too much tension applied to the tape, or from misplacement of the sling or placing the sling too tightly.

- Pain with intercourse (dyspareunia) and loss of sexual function (apareunia), which may be ongoing and may not resolve in some patients.
- Excessive contraction or shrinkage of the tissue surrounding the mesh, and vaginal scarring from causes which include, but are not limited to, chronic inflammation, mesh exposure.
- Neuromuscular problems, including acute and/or chronic pain in the groin, thigh, leg, pelvic and/or abdominal area, and leg weakness, may occur.
- Recurrence of incontinence.
- Bleeding including hemorrhage, or hematoma.
- One or more revision surgeries may be necessary to treat these adverse reactions. Revision surgeries may not resolve complications and are associated with a risk of adverse reactions.
- PROLENE Mesh is a permanent implant that integrates into the tissue. In cases in which the PROLENE Mesh needs to be removed in part or whole, significant dissection may be required. Removal of the implant in whole or in part will not necessarily alleviate the patient's symptoms. Removal of part of the implant can be difficult. Surgery to remove the whole or part of an implant can result in further scarring and tissue damage which, in turn, may have adverse outcomes including severe chronic pain which may not be able to be satisfactorily treated. Surgery to remove the whole or part of the implant may also result in recurrence of SUI. Removal of the eroded mesh will not necessarily prevent further erosions or other adverse events.

Other Adverse Reactions / Undesirable Side Effects

- Seroma
- Urge incontinence, including de novo urge incontinence
- Urinary frequency
- Urinary retention
- Adhesion formation
- Atypical vaginal discharge
- Exposed mesh may cause pain or discomfort to the patient's partner during intercourse.
- Death

The surgeon should convey adverse reactions, undesirable side effects and risks associated with the product and the procedure to the patient prior to undertaking SUI surgery and advise the patient to contact a surgeon in case of any deviation from the normal post-operative course. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

Magnetic Resonance Imaging (MRI)

GYNECARE TVT ABBREVO is Magnetic Resonance (MR) Safe.

Performance(s) / Action(s)

Animal studies show that implantation of PROLENE Mesh elicits an inflammatory reaction in tissues and stimulates the deposition of a thin fibrous layer of tissue that can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The material is not absorbed, nor is it subject to degradation or weakening by the action of tissue enzymes. The device is non-absorbable; therefore, it remains for the lifetime of the patient unless removed.

Sterility

The GYNECARE TVT ABBREVO Continence System is sterile if kept in original, unopened packaging. DO NOT RESTERILIZE. DO NOT REUSE. Do not use if package is opened or damaged.

How Supplied

GYNECARE TVT ABBREVO Continence System is provided sterile (ethylene oxide) for single use.

Storage

Store at or below 25°C (77°F). Do not use after expiry date.

Traceability

A traceability label that identifies the type, size, expiry date, and batch code of the device is included in every package. This label should be affixed or electronically added to the patient's permanent medical record to clearly identify the device that was implanted.

Symbols used on labeling

	Do not re-use		Do not use if package is damaged
	Do not re-sterilize		Sterilized using ethylene oxide
	Caution		Batch code
	Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.		Manufacturer
	Use-by date		Upper limit of temperature
			MR Safe

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