Gynecare

TVT EXACT™

Continence System

GYNECARE TVT EXACT™ kontinenssystem

GYNECARE TVT EXACT™-continentiesysteem

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












9

10

Please read all information carefully.

Failure to properly follow instructions may result in improper functioning of the devices and/or may lead to injury. Properly following instructions will not eliminate the risk of adverse events.

Introduction

This package insert is designed to provide instructions for use of GYNECARE TVT EXACT™ Continen System and Reusable Rigid Catheter Guide. It is not a comprehensive reference to surgical technique for correcting Stress Urinary Incontinence (SUI). This is not a technique manual nor a substitute for appropriate training and experience in surgical technique for correcting SUI. The device is intended to be used only by surgeons trained in the surgical treatment of SUI and specifically in the use of the GYNECARE TVT EXACT Continen System. Training should be in line with the training guidelines specified by the International Urogynecological Association (IUGA) or by other international, national, or local Institutions. Adverse reactions may occur notwithstanding the training and experience of the surgeon. These instructions are recommended for general use of the GYNECARE TVT EXACT Continen System. The removal of the product from the patients (if needed) should be performed by a trained urogynecologist, urologist or gynecologist experienced in the removal of pelvic floor meshes. Variations in use may occur in specific procedures due to individual technique and patient anatomy.

Description

The GYNECARE TVT EXACT Continen System is a sterile, single patient-use procedure kit consisting of:

A. GYNECARE TVT EXACT Continen System Trocar Sheath / Implant Assembly (see Figure 1):

1. Implant
2. Implant Sheath
3. Trocar Sheath
4. Trocar Sheath Cut-out

The GYNECARE TVT EXACT Continen System Trocar Sheath / Implant Assembly consists of one piece of blue (up to 0.28 weight percent phthalocyanine blue, Color Index number 74160) PROLENE™ Polypropylene Mesh (Implant) approximately ½ inch x 18 inches (1.1 cm x 45 cm), covered by clear plastic Sheaths and overlapping in the middle, and held between two white Trocar Sheaths, which are bonded to the Implant and Implant Sheath. Based upon the maximum implanted mesh length, the maximum amount of dye (phthalocyanine blue, Color Index number 74160) is 0.00077 g. PROLENE Polypropylene Mesh is constructed of knitted filaments of extruded polypropylene strands identical in composition to that used in PROLENE™ Polypropylene Nonabsorbable Surgical Sutures. The Implant is approximately 0.027 inches (0.7 mm) thick. This material, when used as a suture, has been reported to be non-reactive and to retain its strength indefinitely in clinical use. PROLENE Mesh is knitted by a process which interlinks each fiber junction.

B. GYNECARE TVT EXACT Continen System Trocar (see Figure 2):

5. Trocar Handle
6. Trocar Sheath Lock
7. Trocar Shaft

The GYNECARE TVT EXACT Continen System Trocar consists of the stainless steel Trocar Shaft and the plastic Trocar Handle. The Trocar Shaft is designed to fit inside the white Trocar Sheaths on the GYNECARE TVT EXACT Continen System Implant / Trocar Sheath Assembly and is used to position the GYNECARE TVT EXACT Continen System Implant in the patient from a vaginal incision up through the abdominal wall.

GYNECARE TVT Reusable Rigid Catheter Guide

(Available separately – not included in GYNECARE TVT EXACT Continen System)

The GYNECARE TVT Rigid Catheter Guide is a non-sterile reusable stainless steel instrument intended to facilitate the identification of the urethra and the bladder neck during the surgical procedure. It is inserted into a Foley catheter (recommended size 18 French) positioned in the bladder via the urethra. To facilitate insertion, it can be lubricated with gel. The Rigid Catheter Guide must be cleaned and sterilized prior to use. Instructions can be found in the Instructions for Cleaning Reusable Instrument section.

Indications / Intended Use

The GYNECARE TVT EXACT Continen System is intended to be used as a pubo-urethral sling for treatment of female Stress Urinary Incontinence, resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

The GYNECARE TVT EXACT Continen System Trocar is a single use device intended to aid in the placement of the GYNECARE TVT EXACT Continen System retropubically.

The GYNECARE TVT Rigid Catheter Guide is a reusable device intended to facilitate the identification of the urethra and bladder neck during the placement of the GYNECARE TVT EXACT Continen System.

Contraindications

As with any suspension surgery, this procedure should not be performed in pregnant patients. Additionally, because the PROLENE Polypropylene Mesh will not stretch significantly, this procedure should not be performed in patients with future growth potential, including women with plans for future pregnancy.

Patient Target Group and Factors

GYNECARE TVT EXACT is intended for use for non-pregnant adult female patients not planning future pregnancies, who are affected by SUI resulting from urethral hypermobility and/or intrinsic sphincter deficiency. Surgeons should use their surgical experience and judgment to determine if PROLENE Mesh is appropriate for certain patients. Patient-specific factors may impair wound healing, which may increase the likelihood of adverse reactions.

Application / Instructions for Use

Note: Hand positions used in surgery may vary from those shown in the illustrations.

1. The procedure can be carried out under local anesthesia, but it can also be performed using regional or general anesthesia.
2. Before the patient is prepped and draped, she should be placed in the lithotomy position, taking care to avoid hip flexion greater than 60 degrees.
3. Insert an 18 French Foley catheter into the bladder and leave it to open drainage.
4. At the level of the mid-urethra, inject a small amount of local anesthesia submucosally to create a space between the vaginal wall and the peri-urethral fascia. The extent of dissection required for placement is minimal. Only a small para-urethral incision is required over the mid-urethra to position the tip of the Trocar Sheath. Using forceps, grasp the vaginal wall at each side of the urethra. Using a small scalpel, make a sagittal incision no more than 1.5 cm long starting approximately 1.0 cm cephalad from the urethral meatus. This incision will be positioned over the mid-urethral zone and will allow for subsequent passage of the Implant.
5. With a small pair of blunt scissors, make two small para-urethral lateral dissections (approximately 0.5 cm to 1.0 cm) to accommodate the tips of the Trocar Sheaths.
6. Identify the two Trocar Sheath exit sites, which should be 2 cm to 2.5 cm on each side of the midline, immediately above the pubic symphysis (See Figure 3). Either mark these sites or, if desired, place two small 3 mm to 4 mm transverse stab incisions at the intended exit site. To avoid the inferior epigastric vessels, it is important that the intended exit sites be not more than 2.5 cm from the midline. It is important that the exit sites for the Trocar Sheath passages be near the midline and close to the superior aspect of the pubic bone in order to avoid anatomic structures in the abdomen, inguinal area and lateral pelvic sidewall.
7. If retroperitoneal infiltration of local anesthesia is not performed then consider infiltrating the retropubic space with two injections of normal saline on either side of midline. Starting at the Trocar exit sites, pass an 18-gauge needle along the back of the pubic bone until the tip of the Trocar touches the endopelvic fascia. As the Trocar is withdrawn, inject 30 mL to 50 mL. By so doing, it opens up the retropubic space to further minimize the risk of bladder puncture during retropubic Trocar passage.
8. Confirm that all urine has been drained from the bladder. Once the bladder is empty, insert the GYNECARE TVT Reusable Rigid Catheter Guide (available separately) into the channel of the 18 French Foley catheter. The handle of the GYNECARE TVT Reusable Rigid Catheter Guide should be fixed around the catheter at its proximal end. The purpose of placing the GYNECARE TVT Reusable Rigid Catheter Guide into the catheter is to allow contralateral displacement of the bladder, bladder neck and urethra away from the tip of the Trocar Sheath as it is passes through the retropubic space.

9. Place the Trocar Shaft inside one of the two white Trocar Sheaths (See Figure 4). Secure the Trocar Sheath to the Trocar Handle by hooking the Trocar Sheath Cut-out onto the Trocar Sheath Lock on the Trocar Handle (See Figure 5). NOTE: Ensure that the Trocar Sheath Cut-out goes completely over the Trocar Sheath Lock and is holding the Trocar Sheath on the Trocar Shaft securely. Be careful not to manipulate the Trocar Sheath appendage hanging past the Trocar Sheath Lock during the procedure, as that may result in the unintended disengagement of the Trocar Sheath Lock.
10. Using the Rigid Catheter Guide, gently push the tip of the 18 French Foley catheter toward the posterior lateral wall of the bladder opposite to the intended Trocar Sheath passage (see Figure 6). For example, by pushing toward the patient's left side the bladder will go from a spherical to a spheroid configuration. This moves the bladder away from the back of the pubic symphysis (see Figure 7A and 7B). Additionally, it moves the bladder neck and the urethra to the left as the Trocar Sheath is passed on the patient's right side, thereby minimizing the risk of bladder perforation (see Figure 8A and 8B).
11. Hold the Trocar Handle using your dominant hand. Pass the tip of the white Trocar Sheath that has been mounted on the Trocar Shaft (see previous Step 8), para-urethrally through the urogenital diaphragm at the level of the mid-urethra. Initial insertion of the device is controlled by using the tip of the index finger of the non-dominant hand, which is placed in the vagina under the anterior vaginal wall, just lateral to the sub-urethral incision. The curved part of the Trocar Shaft should rest in the palm of the non-dominant hand (See Figure 9). Pass the Trocar Sheath through the urogenital diaphragm into the retropubic space. During the initial placement into the para-urethral dissected space, the Trocar Sheath tip should be oriented horizontally, i.e., in the frontal plane. During passage through the urogenital diaphragm, lower the Trocar Handle to ensure that the Trocar Sheath Tip passes vertically while staying in close contact to the back of the pubic symphysis. After passage through the urogenital diaphragm, resistance to the passage of the Trocar Sheath is significantly reduced once it enters the retropubic space.
12. At this point, the non-dominant hand is moved from the vagina to the suprapubic exit site. The Trocar Sheath tip is guided through the retropubic space staying as close to the back of the pubic symphysis as possible. This is achieved by lowering the Trocar Handle, thereby pressing the Trocar Sheath tip against the back of the pubic bone.
13. During passage through the retropubic space, aim the Trocar Sheath tip toward the pre-marked abdominal exit site.
14. Move the Trocar Sheath tip upward toward the abdominal skin exit sites, keeping in close contact with the pubic bone until exiting the skin (See Figure 10). Once the Trocar Sheath tip exits the skin, grasp the exposed Trocar Sheath tip with a clamp. Release the Trocar Sheath from the Trocar Sheath Lock on the Trocar Handle by pushing the Trocar Sheath appendage laterally and off the Trocar Sheath Lock, and carefully withdraw the Trocar Shaft from within the Trocar Sheath. DO NOT PULL the Trocar Sheath up any further.
15. The procedure is now repeated on the patient's other side while repeating steps 9 to 14. NOTE: IN ORDER TO MINIMIZE THE RISK OF BLADDER INJURY, IT IS IMPORTANT THAT THE BLADDER NOW BE DISPLACED TO THE CONTRALATERAL SIDE USING THE MANEUVERS OUTLINED IN STEP 10.
16. Once both white Trocar Sheaths have been passed and before the Implant is pulled into place, remove the 18 French Foley catheter and perform a cystoscopy (70 degree lens recommended).
17. Once bladder integrity is confirmed, gently pull the Trocar Sheaths upward to bring the Implant loosely (i.e., without tension) under the mid-urethra. Cut the Implant bilaterally close to the connection to the Trocar Sheaths. Adjust the Implant so that leakage is reduced, allowing only a few drops of urinary leakage to occur under stress. For this, use patient feedback, i.e., coughing with a full bladder (approximately 300 mL). In case of general anesthesia, verify tension of the tape exerting manual pressure on the suprapubic area with a full bladder (approximately 300 mL).
18. Grasp the Implant Sheaths that surround the Implant with clamps, taking care not to grasp the Implant. Then individually remove the Implant Sheaths by gently pulling up on the clamps, away from the abdomen, one at a time. To avoid putting tension on the Implant, a blunt instrument (scissors or forceps) should be placed between the urethra and the Implant during removal of the Implant Sheaths.

NOTE: Premature removal of the Sheaths may make subsequent adjustments difficult.

19. After proper adjustment of the Implant, close the vaginal incision. The abdominal ends of the Implant are then cut and left in the subcutis. Do not suture the Implant.
20. Close the skin incisions with suture or surgical skin adhesive.
21. Empty the bladder. Following this procedure, postoperative catheterization is not typically required. The patient should be encouraged to try to empty their bladder 2 to 3 hours after the operation.

Disposal

Discard unintentionally opened / partially used / used devices and packaging according to your facility's policies and standard operating procedures concerning biohazardous materials and waste.

Warnings and Precautions

- Do not use GYNECARE TVT EXACT Continence System in patients who are on anti-coagulation therapy.
- Do not use the GYNECARE TVT EXACT Continence System in a patient who has a urinary tract infection.
- Surgeons should be familiar with surgical technique for mid-urethral sling placement and should be adequately trained in the GYNECARE TVT EXACT Continence System implantation procedure before employing the GYNECARE TVT EXACT Continence System. It is important that the Implant be located without tension under mid-urethra.
- Acceptable surgical practice should be followed for the procedure as well as for the management of contaminated or infected wounds.
- The procedure should be performed with care to avoid large vessels, nerves, bladder, and bowel. Attention to local anatomy and proper passage of the Trocar Sheaths may minimize certain risks.
- Retropubic bleeding may occur post-operatively. Observe for any symptoms or signs before releasing the patient from the hospital.
- Cystoscopy should be performed to confirm bladder integrity or recognize a bladder perforation.
- The Rigid Catheter Guide should be gently pushed into the Foley catheter so that the Catheter Guide does not extend into the holes of the Foley catheter.
- When removing the Rigid Catheter Guide, open the handle completely so that the Foley catheter remains properly in place.
- Do not remove the plastic Sheath until the Implant has been properly positioned.
- Ensure that the Implant is placed with minimal tension under the mid-urethra.
- Do not perform this procedure if you think the surgical site may be infected or contaminated.
- Since limited clinical information is available about pregnancy following pubo-urethral sling procedure with the GYNECARE TVT EXACT Continence System, the patient should be counseled that future pregnancies may negate the effects of the surgical procedure and the patient may again become incontinent.
- Since limited clinical information is available with delivery following a pubo-urethral procedure with the GYNECARE TVT EXACT Continence System, in case of pregnancy, mode of delivery should be determined by the obstetrician in consultation with the surgeon.
- Post-operatively, the patient is recommended to refrain from heavy lifting and/or exercise (e.g., cycling, jogging) for at least four to six weeks and intercourse for at least one month. The patient can return to other normal activity after one or two weeks.
- Should dysuria, bleeding, or other problems occur, the patient should be instructed to contact the surgeon immediately.
- All surgical instruments are subject to wear and damage under normal use. Before use, the instrument should be visually inspected. Defective instruments or instruments that appear to be corroded should not be used and should be discarded.
- As with other incontinence procedures, de novo detrusor instability may occur following a pubo-urethral procedure utilizing the GYNECARE TVT EXACT Continence System. To minimize this risk, make sure to place the Implant tension-free in the mid-urethral position, as described above.
- Do not contact the PROLENE Mesh with any staples, clips or clamps, as mechanical damage to the mesh may occur.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and surgeons.
- Prophylactic antibiotics can be administered according to the surgeon's usual practice.
- In the event of a damaged product, the product should be discarded and a new one retrieved to complete the procedure. In the event of a product malfunction during use, continue or discontinue use at surgeon's discretion.

Adverse Reactions / Undesirable Side Effects

- Punctures or lacerations of vessels, nerves, structures, or organs, including the bladder, urethra, or bowel, may occur and may require surgical repair.
- Local irritation at the wound site may occur.
- As with any implant, a foreign body response will occur, the extent of which may differ. This response could result in extrusion, erosion, exposure, fistula formation and/or chronic inflammation, the severity of which is unpredictable, or other adverse reactions, which may be ongoing.
- Fistula formation, acute and chronic inflammation and ongoing risk of mesh extrusion, exposure, or erosion into the vagina or other structures or organs (such as bladder, urethra or rectum), which may be difficult to treat and result in consequent pain. Mesh extrusion, exposure, or erosion into the vagina may also cause offensive vaginal discharge.

- Infection following transvaginal implantation. As with all surgical procedures and the implantation of foreign bodies, there is a risk of infection and PROLENE Mesh may potentiate an existing infection.
- Pain – which may be severe and chronic
- Temporary or chronic voiding dysfunction (or difficulty voiding) or urinary retention/obstruction independent from that cause by overcorrection or urethral hypermobility, i.e., too much tension applied to the tape, or from misplacement of the sling or placing the sling too tightly.
- Pain with intercourse (dyspareunia), and loss of sexual function (apareunia), which may be ongoing and may not resolve in some patients.
- Excessive contraction or shrinkage of the tissue surrounding the mesh, and vaginal scarring from causes including, but not limited to, chronic inflammation, mesh exposure.
- Neuromuscular problems, including acute and/or chronic pain in the groin, thigh, leg, pelvic and/or abdominal area, and leg weakness, may occur.
- Recurrence of incontinence.
- Bleeding, including hemorrhage, or hematoma.
- One or more revision surgeries may be necessary to treat these adverse reactions. Revision surgeries may not resolve complications and are associated with a risk of adverse reactions.
- PROLENE Mesh is a permanent implant that integrates into the tissue. In cases in which the PROLENE Mesh needs to be removed in part or whole, significant dissection may be required. Removal of the implant in whole or in part will not necessarily alleviate the patient's symptoms. Removal of part of the implant can be difficult. Surgery to remove the whole or part of an implant can result in further scarring and tissue damage which, in turn, may have adverse outcomes including severe chronic pain which may not be able to be satisfactorily treated. Surgery to remove the whole or part of the implant may also result in recurrence of SUI. Removal of the eroded mesh will not necessarily prevent further erosions or other adverse events.

Other Adverse Reactions / Undesirable Side Effects

- Seroma
- Urge incontinence, including de novo urge incontinence
- Urinary frequency
- Urinary retention
- Adhesion formation
- Atypical vaginal discharge
- Exposed mesh may cause pain or discomfort to the patient's partner during intercourse.
- Death

The surgeon should convey adverse reactions, undesirable side effects and risks associated with the product and the procedure to the patient prior to undertaking SUI surgery and advise the patient to contact a surgeon in case of any deviation from the normal post-operative course. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

Magnetic Resonance Imaging (MRI)

GYNECARE TVT EXACT Continence System is Magnetic Resonance (MR) Safe.

Performance(s) / Action(s)

Animal studies show that implantation of PROLENE Mesh elicits an inflammatory reaction in tissues and stimulates the deposition of a thin fibrous layer of tissue that can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The material is not absorbed, nor is it subject to degradation or weakening by the action of tissue enzymes. The device is non-absorbable; therefore, it remains for the lifetime of the patient unless removed.

Sterility

The GYNECARE TVT EXACT Continence System is sterile if kept in original, unopened packaging. DO NOT RESTERILIZE. DO NOT REUSE. Do not use if package is opened or damaged.

The reusable GYNECARE TVT Rigid Catheter Guide is supplied separately and is non-sterile. It is to be cleaned and sterilized prior to each use as described below.

Instructions for Cleaning Reusable Instrument GYNECARE TVT Rigid Catheter Guide (available separately)

To ensure the reliability and functionality of the GYNECARE TVT Rigid Catheter Guide, clean and sterilize the instrument before initial use and after each procedure. The following are suggested pre-cleaning and automated methods for cleaning the instrument. Each instrument must be treated individually, and matt surfaces must be protected against scratches. The design of this instrument permits effective cleaning and reliable sterilization. The user must ensure that cleaning and sterilization are performed in accordance with guidelines, standards and National Health Authority requirements.

Point-of-Use Care of the Instrument

Wipe the instrument clean of visible soil. Properly contain the instrument for transport to a dedicated cleaning work area, ensuring the instrument is kept moist to avoid drying of soil.

Pre-Cleaning of the Instrument

1. The section of the Rigid Catheter Guide that clips onto the long shaft should be opened so that all surfaces are easily accessible for cleaning.
2. Rinse the instrument thoroughly with tap water.
3. Prepare a neutral pH enzymatic solution suitable for stainless steel instruments per the manufacturer's recommendations.
4. Soak the instrument components in the enzymatic cleaner for a minimum of 10 minutes.
5. Wash in the enzymatic solution at a temperature recommended by the manufacturer. Remove any contamination from body fluids or tissues. The Rigid Catheter Guide should be cleaned using a soft-bristled brush (Spectrum M16 or equivalent) and scrubbing the surfaces for a minimum of 2 minutes.
6. Remove the instrument from the enzymatic solution and rinse thoroughly under running tap water for a minimum of 1 minute.
7. Place the instrument in an ultrasonic cleaner with fresh enzymatic detergent for a minimum of 10 minutes. The frequency should be 25 kHz to 50 kHz (target 35 kHz).
8. Remove the instrument from the ultrasonic cleaner and rinse thoroughly under running tap water for a minimum of 1 minute.
9. Visually inspect the instrument for cleanliness. If necessary, repeat steps 1 to 8. Following pre-cleaning, the instrument should be cleaned in a washer-disinfector that conforms to ISO 15883-1 and 2.

Automated Cleaning Method

After the pre-cleaning steps have been carried out, the instruments can be cleaned in an automated washing cycle. Automatic washing cycles are suitable for stainless steel instruments. One recommended cycle is described below:

PHASE	MINIMUM TIME (MINUTES)	TEMPERATURE AND WATER QUALITY	DETERGENT TYPE AND CONCENTRATION
Pre-wash	02:00	Cold tap water	N/A
Enzyme wash	15:00	Hot tap water or detergent manufacturer's recommendation Minimum 176°F (80°C)	Enzymatic Detergent (Enzol 2 oz./gallon or an equivalent product)
Rinse 1	02:00	Hot tap water	N/A
Rinse 2	15:00	Hot tap water	N/A
Final rinse	05:00	Hot tap water	N/A
Rinse with demineralized or purified water	02:00	Demineralized or purified water 176°F (80°C)	N/A
Thermal disinfection (optional)	01:00	Demineralized or purified water 194°F (90°C)	N/A
Drying	10:00-15:00	199.4°F (93°C)	N/A

Thermal Disinfection (Optional)

Note that disinfection is mandatory in some areas of the world, e.g., EU. The user must ensure that cleaning, disinfection, and sterilization are conducted in accordance with the appropriate guidelines, standards, or National Health Authority requirements.

Thermal disinfection can be included as part of the automatic cleaning cycle. The instrument can be thermally disinfected at 194°F (90°C) or greater for a minimum of 1 minute to render the instrument safe for handling.

Note: Do not exceed 212°F (100°C) for longer than 10 minutes.

Sterilization Recommendations for Reusable Instrument GYNECARE TVT Rigid Catheter Guide (available separately)

The GYNECARE TVT Rigid Catheter Guide is supplied non-sterile. To sterilize, steam autoclave prior to each use. Steam (moist heat) sterilization is recommended in an approved, pre-vacuum cycle listed below. The steam sterilizer must be validated to the requirements of any local standards and guidance such as EN285 or AAMI/ANSI ST8. Ensure that the steam sterilizer cycle is chosen that is designed to remove air from porous loads and/or lumened instruments in accordance to manufacturer's instructions.

Prior to steam sterilization, ensure the instrument is dry and package the instrument in two layers of sterilization wrap or other suitable sterile barrier package. The sterile barrier packaging should meet the requirements in ISO 11607-1.

WARNING: The section of the Rigid Catheter guide that clips onto the long shaft must be open during cleaning and sterilization. If it is not open during cleaning and sterilization, the cleaning process may be less effective, and sterility of the instrument may be compromised.

The types of sterilization equipment and available cycles may vary by geographic region. The following sterilization parameters have been validated to deliver a sterility assurance level of 10⁻⁶ for this instrument:

Market Area	Cycle Type	Temperature Set Point	Exposure Time (minutes)	Minimum Drying Time (minutes)
USA and its territories	Pre-vacuum	270°F/132°C	4	35
EU and others	Pre-vacuum	273.2°F/134°C	3	35

Note:

The above table includes the minimum temperature and cycle time validated to assure sterility.

It is the responsibility of the end user to assure sterility of the product when using sterilization process recommended, since bioburden and sterilization equipment will vary.

GYNECARE TVT Rigid Catheter Guide Maintenance

Before each use, inspect the instrument. Check to ensure that the long end which traverses the catheter channel has no sharp edges or burrs.

How Supplied

The GYNECARE TVT EXACT Continence System is provided sterile (ethylene oxide) for single use. The reusable GYNECARE TVT Rigid Catheter Guide is supplied separately and is non-sterile. The accessory is to be cleaned and sterilized prior to each use as described above.

Storage

Store at or below 25°C (77°F). Do not use after expiry date.

Traceability

A traceability label that identifies the type, size, expiry date, and batch code of the device is included in every package. This label should be affixed or electronically added to the patient's permanent medical record to clearly identify the device that was implanted.

Symbols Used on Labeling

 Caution	 Do not use if package is damaged	 Manufacturer
 Do not re-use	 Use-by date	 Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
 Do not re-sterilize	 Sterilized using ethylene oxide	 Upper limit of temperature
 Catalogue number	 Batch code	 MR Safe

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© Ethicon, Inc. 2021