TVT™ Obturator System

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

CE 2797


Made in Switzerland
© Ethicon, Inc. 2021
P18070
LAB0010875v7
01/2021


Ethicon, SÄRL
Puits Godet 20
Neuchâtel
CH-2000
Switzerland
1-877-ETHICON
+1-513-337-6928


1


2


3


4


5


6


7


8

9

10

11

GYNECARE TVT™ Obturator System Tension-free Support for Incontinence

GYNECARE TVT™ Obturator Device, Sterile Single Use

GYNECARE TVT™ Obturator Helical Passers, Sterile Single Use

GYNECARE TVT™ Obturator Atraumatic Winged Guide, Sterile Single Use

Please read all information carefully.

Failure to properly follow instructions may result in improper functioning of the device and may lead to injury. Properly following instructions will not eliminate the risk of adverse events.

Introduction

This package insert is designed to provide instructions for use of the GYNECARE TVT Obturator System, including the GYNECARE TVT Obturator Device, Helical Passers and Atraumatic Winged Guide. It is not a comprehensive reference to surgical technique for correcting SUI (Stress Urinary Incontinence). This is not a technique manual nor a substitute for appropriate training and experience in surgical technique for correcting SUI. The device is intended to be used only by surgeons trained in the surgical treatment of SUI and specifically in implanting the GYNECARE TVT Obturator Device. Training should be in line with the training guidelines specified by the International Urogynecological Association (IUGA) or by other international, national, or local Institutions. Adverse reactions may occur notwithstanding the training and experience of the surgeon. These instructions are intended for general use of the device. The removal of the product from the patients (if needed) should be performed by a trained urogynecologist, urologist or gynecologist experienced in the removal of pelvic floor meshes. Variations in use may occur in specific procedures due to individual technique and patient anatomy.

Description

The GYNECARE TVT Obturator System is a sterile, single patient use procedure kit consisting of:

GYNECARE TVT Obturator Device

The GYNECARE TVT Obturator Device consists of one piece of undyed or blue (up to 0.28 weight percent phthalocyanine blue, Color Index number 74160) PROLENE™ polypropylene mesh (tape) approximately 1/2 inch x 18 inches (1.1 cm x 45 cm) covered by clear plastic sheath overlapping in the middle. Plastic tube receptacles are attached at each end. Based upon the maximum implanted mesh length, the maximum amount of dye (phthalocyanine blue, Color Index number 74160) is 0.00067 g. The GYNECARE TVT Obturator Device is available in either mechanical cut or laser cut versions for the surgeon's preference. To determine if the GYNECARE TVT Obturator System implant is mechanical or laser cut, consult the product code on the device packaging; an (L) at the end of the number indicates the laser cut mesh. PROLENE Polypropylene Mesh is constructed of knitted filaments of extruded polypropylene strands identical in composition to that used in PROLENE Polypropylene non-absorbable surgical suture. The Implant is approximately 0.027 inches (0.7 mm) thick. This material, when used as a suture, has been reported to be non-reactive and to retain its strength indefinitely in clinical use. PROLENE Mesh is knitted by a process that interlinks each fiber junction.

GYNECARE TVT Helical Passers

The GYNECARE TVT Helical Passers are two stainless steel, curved wire passers with plastic handles that are designed to deliver the GYNECARE TVT Obturator Device. Helical Passers are provided as left and right units, pre-assembled to the GYNECARE TVT Obturator Device. The Helical Passer MUST not be bent or deformed in any way.

GYNECARE TVT Atraumatic Winged Guide

The GYNECARE TVT Atraumatic Winged Guide is a stainless steel accessory instrument, which facilitates the passage of the GYNECARE TVT Helical Passers through the dissection tract.

Indications / Intended Use

The GYNECARE TVT Obturator Device is intended to be used in women as a sub-urethral sling for the treatment of stress urinary incontinence (SUI) resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

The GYNECARE TVT Obturator Helical Passers and Atraumatic Winged Guide are intended to aid in the placement of the GYNECARE TVT Obturator Device.

Contraindications

As with any suspension surgery, this procedure should not be performed in pregnant patients. Additionally, because the PROLENE Polypropylene mesh will not stretch significantly, this procedure should not be performed in patients with future growth potential, including women with plans for future pregnancy.

Patient Target Group and Factors

GYNECARE TVT Obturator is intended for use for non-pregnant adult female patients not planning future pregnancies, who are affected by SUI resulting from urethral hypermobility and/or intrinsic sphincter deficiency.

Surgeons should use their surgical experience and judgment to determine if PROLENE Mesh is appropriate for certain patients. Patient-specific factors may impair wound healing, which may increase the likelihood of adverse reactions.

Application / Instructions for Use

Note: Hand positions used in surgery may vary from those shown in illustrations.

1. Place the patient in the lithotomy position with the hips hyperflexed over the abdomen. The tip of the coccyx should be positioned flush with the edge of the table. Abduction of the legs should be performed after flexion of the hips is realized.
2. The procedure can be carried out under local, regional, or general anesthesia.
3. If desired, retract the labia to provide additional exposure.
4. Insert a urethral catheter into the bladder and empty the bladder.
5. Mark the exit points of the plastic tubes by tracing a horizontal line at the level of the urethral meatus, and a second line parallel and 2 cm above the first line. Locate the exit points on this line, 2 cm lateral to the folds of the thigh (the skin may be flattened by stretching). Mark the exit points, alternatively a 5 to 10 mm incision may be made at each exit point or at a later stage of the procedure. *(See Figure 1)*
6. Using Allis clamps for traction, make a 1 cm midline incision in the vaginal mucosa starting 1 cm proximal to the urethral meatus and two lateral sharp dissections.

(Note: It is suggested that the device insertion be completed on one side before beginning dissection of the second side.)

After initiating sharp dissection, continue by using a “push-spread technique” to perform blunt dissection, preferably using pointed, curved scissors. The path of the lateral dissection should be oriented at a 45° angle from the midline, with the scissors oriented on the horizontal plane. *(See Figure 2)* Continue dissection towards the junction between the body of the pubic bones and the inferior pubic ramus. *(See Figure 2)*

- 2A. Horizontal Plane
- 2B. Ischio-Pubic Ramus
- 2C. Midline Plane
- 2D. Vaginal Dissection

When the junction between the body of the pubic bones and the inferior pubic ramus is reached, perforate the obturator membrane. A loss of resistance can be felt when the membrane is perforated. The channel should be approximately 5 to 7 mm in diameter and no deeper than 5 cm. Dissection beyond 5 cm may allow unintended entry into the Space of Retzius. If the bone is not reached after dissecting 5 cm, re-evaluate that the angle of dissection is correct.

7. Remove the internal package workstation from the external package. Then remove the GYNECARE TVT Winged Guide from the package workstation. *(See Figure 3)*
 - 3A. Winged Guide
 - 3B. Helical Passer/Device Assembly
8. Insert the GYNECARE TVT Winged Guide into the dissected tract until it passes the inferior pubic ramus and enters the opening previously made in the obturator membrane. Loss of resistance can be felt as the Winged Guide passes through the obturator membrane.

If difficulty is encountered during insertion of the guide, reconfirm the direction of the tract with the scissors.

(Note: The open side of the guide must be facing the surgeon. The bendable tab can be bent to increase the length of the guide if needed, see Figure 5.)

9. Remove the GYNECARE TVT Helical Passers/Device Assembly and the GYNECARE TVT Obturator Device assembly from the sterile pack. *(see Figure 3 for components)*

(Note: To ensure correct orientation of the Helical Passers and tape, verify that the GYNECARE logo and thumb indent on the plastic handle are facing the surgeon, and that the points are on the outside facing the surgeon. The Helical Passer in the surgeon's left hand must be used on the patient's right side; see Figure 4.)

10. Place one of the Helical Passers on the sterile drape or other suitable sterile location until needed. Assume that the tape is not twisted.
11. Insert the correct GYNECARE TVT Helical Passer into the dissected tract following the channel of the GYNECARE TVT Winged Guide. Push the device inward, traversing, and slightly passing the obturator membrane. Make sure the device handle is oriented so the straight tip of the Helical Passer is aligned with the channel in the GYNECARE TVT Winged Guide and remains in that orientation until the tip traverses the obturator membrane. *(see Figure 5)*
12. Once in this position, remove the GYNECARE TVT Winged Guide and keep sterile for later use on the same patient.
13. Once the GYNECARE TVT Winged Guide has been removed, rotate the handle of the Helical Passer as you simultaneously move towards the midline until the handle is vertical to the floor. *(see Figure 6)*

- 6A. Starting position
- 6B. Ending position

(Note: Never allow the handle to be oriented horizontal to the floor.)

14. The point of the Helical Passer should exit near the previously determined exit points. *(see Figure 7)* However, slight skin manipulation may be required. If the skin incision has not been previously made, make it at the point where the tip of the Helical Passer tents the skin. When the tip of the plastic tube appears at the skin opening, grasp the pointed tip of the plastic tube with a clamp and, while stabilizing the tube near the urethra with the thumb, remove the Helical Passer by a reverse rotation of the handle. *(see Figure 8)*

15. Pull the plastic tube completely through the skin until the tape appears. (see Figure 9)
16. Repeat the technique on the patient's other side ensuring that the tape lies flat under the urethra. (see Figure 10)

(Note: If a twist in the tape is discovered, ensure that the twist is not positioned under the urethra after the excess tape is pulled through.)

17. When both plastic tubes have been extracted through the skin incisions, cut the plastic tubes from the tape and plastic sheaths. Position the tape loosely, e.g., without tension, and flat under the mid-urethra. At this stage, it is advisable to perform a cough test or Credé maneuver (i.e. exerting manual pressure on the suprapubic area, in particular if the patient is under general anesthesia) after instilling the bladder with 300 mL (or up to the patient's preoperatively determined bladder capacity) to adjust mesh implant tensioning to the individual patient's characteristics. (see Figure 11)

When the tape is in position, remove the plastic sheath that covers the tapes. Place a blunt instrument (e.g., scissors or forceps) between the urethra and the tape during removal of the plastic sheaths, or use other suitable means during sheath removal, to avoid positioning the tape with tension.

(Note: Premature removal of the sheath may make subsequent adjustments difficult.)

18. Following tape adjustment close the vaginal incision. Cut the tape ends at the exit points just below the skin of the inner thigh. Close the skin incisions with suture or surgical skin adhesive.
19. Cystoscopy can be performed at the discretion of the surgeon. If cystoscopy was performed following the first passage, make sure the bladder is emptied prior to initiating passage of the second side. Post-operative indwelling catheterization is not typically required. The patient should be encouraged to try to empty the bladder 2 to 3 hours after the operation.

Disposal

Discard unintentionally opened / partially used / used devices and packaging according to your facility's policies and standard operating procedures concerning biohazardous materials and waste.

Warnings and Precautions

- Do not use GYNECARE TVT Obturator System in patients who are on anti-coagulation therapy.
- Do not use GYNECARE TVT Obturator System in a patient who has a urinary tract infection.
- Surgeons should be familiar with surgical technique for mid-urethral sling and should be adequately trained in the GYNECARE TVT Obturator System procedure before employing the GYNECARE TVT Obturator Device. It is important that the Implant be located without tension under mid-urethra.
- Acceptable surgical practice should be followed for the procedure as well as for the management of contaminated or infected wounds.
- The procedure should be performed with care to avoid large vessels, nerves, bladder, and bowel. Attention to patient anatomy and correct passage of the device may minimize certain risks.
- Bleeding may occur post-operatively. Observe for any symptoms or signs before releasing the patient from hospital.
- Although bladder injury is unlikely to occur with this technique, cystoscopy may be performed at the discretion of the surgeon.
- Do not remove the plastic sheaths until the Implant has been properly positioned.
- Ensure that the tape is placed with no tension under the mid-urethra.
- Do not perform this procedure if you think the surgical site may be infected or contaminated.
- Since no clinical information is available about pregnancy following sub-urethral sling procedure with the GYNECARE TVT Obturator System, the patient should be counselled that future pregnancies may negate the effects of the surgical procedure and the patient may again become incontinent.
- Since limited clinical information is available with delivery following a sub-urethral sling procedure with the GYNECARE TVT Obturator System, in case of pregnancy, mode of delivery should be determined by the obstetrician in consultation with the surgeon.
- Post-operatively, the patient should be advised to refrain from heavy lifting and/or exercise (e.g., cycling, jogging) for at least four to six weeks and intercourse for at least one month. The patient can usually return to other normal activity after one or two weeks.
- Should dysuria, bleeding or other problems occur, the patient should be instructed to contact the surgeon immediately.
- Transient leg pain lasting 24–48 hours may occur and can usually be managed with mild analgesics.
- As with other incontinence procedures, de novo detrusor instability may occur following a sub-urethral sling procedure utilizing the GYNECARE TVT Obturator System. To minimize this risk, make sure to place the tape as described above. Do not contact the PROLENE Mesh with any staples, clips or clamps as mechanical damage to the mesh may occur.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and surgeons.
- Prophylactic antibiotics can be administered according to the surgeon's usual practice.
- In the event of a damaged product, the product should be discarded and a new one retrieved to complete the procedure. In the event of a product malfunction during use, continue or discontinue use at surgeon's discretion.

Adverse Reactions / Undesirable Side Effects

- Punctures or lacerations of vessels, nerves, structures, or organs, including the bladder, urethra, or bowel, may occur and may require surgical repair.
- Local irritation at the wound site may occur.
- As with any implant, a foreign body response will occur, the extent of which may differ. This response could result in extrusion, erosion, exposure, fistula formation and/or chronic inflammation, the severity of which is unpredictable or other adverse reactions, which may be ongoing.
- Fistula formation, acute and chronic inflammation and ongoing risk of mesh extrusion, exposure, or erosion into the vagina or other structures or organs (such as bladder, urethra or rectum) which may be difficult to treat and result in consequent pain. Mesh extrusion, exposure, or erosion into the vagina may also cause offensive vaginal discharge.
- Infection following transvaginal implantation. As with all surgical procedures and the implantation of foreign bodies, there is a risk of infection and PROLENE Mesh may potentiate an existing infection.
- Pain - which may be severe and chronic.
- Temporary or chronic voiding dysfunction (or difficulty voiding) or urinary retention/obstruction independent from that caused by overcorrection or urethral hypermobility, i.e., too much tension applied to the tape, or from misplacement of the sling or placing the sling too tightly.
- Pain with intercourse (dyspareunia), and loss of sexual function (apareunia), which may be ongoing and may not resolve in some patients.
- Excessive contraction or shrinkage of the tissue surrounding the mesh, and vaginal scarring from causes which include, but are not limited to, chronic inflammation, mesh exposure.
- Neuromuscular problems, including acute and/or chronic pain in the groin, thigh, leg, pelvic and/or abdominal area, and leg weakness, may occur.
- Recurrence of incontinence.
- Bleeding, including hemorrhage, or hematoma.
- One or more revision surgeries may be necessary to treat these adverse reactions. Revision surgeries may not resolve complications and are associated with a risk of adverse reactions.
- PROLENE Mesh is a permanent implant that integrates into the tissue. In cases in which the PROLENE Mesh needs to be removed in part or whole, significant dissection may be required. Removal of the implant in whole or in part will not necessarily alleviate the patient's symptoms. Removal of part of the implant can be difficult. Surgery to remove the whole or part of an implant can result in further scarring and tissue damage which, in turn, may have adverse outcomes including severe chronic pain which may not be able to be satisfactorily treated. Surgery to remove the whole or part of the implant may also result in recurrence of SUI. Removal of the eroded mesh will not necessarily prevent further erosions or other adverse events.

Other Adverse Reactions / Undesirable Side Effects

- Seroma
- Urge incontinence, including de novo urge incontinence
- Urinary frequency
- Urinary retention
- Adhesion formation
- Atypical vaginal discharge
- Exposed mesh may cause pain or discomfort to the patient's partner during intercourse.
- Death

The surgeon should convey adverse reactions, undesirable side effects and risks associated with the product and the procedure to the patient prior to undertaking SUI surgery and advise the patient to contact a surgeon in case of any deviation from the normal post-operative course. Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

Magnetic Resonance Imaging (MRI)

GYNECARE TVT Obturator System is Magnetic Resonance (MR) Safe.

Performance(s) / Action(s)

Animal studies show that implantation of PROLENE Mesh elicits an inflammatory reaction in tissues and stimulates the deposition of a thin fibrous layer of tissue that can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The material is not absorbed, nor is it subject to degradation or weakening by the action of tissue enzymes. The device is non-absorbable, therefore it remains for the lifetime of the patient unless removed.

Sterility

The GYNECARE TVT Obturator Continence System is sterile if kept in original, unopened packaging. DO NOT RESTERILIZE. DO NOT REUSE. Do not use if package is opened or damaged.

How Supplied

The GYNECARE TVT Obturator System is provided sterile (ethylene oxide) for single use.

Storage

Store at or below 25°C (77°F). Do not use after expiry date.

Traceability

A traceability label that identifies the type, size, expiry date, and batch code of the device is included in every package. This label should be affixed or electronically added to the patient's permanent medical record to clearly identify the device that was implanted.

Symbols Used On Labeling

	Do not re-use		Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	Do not re-sterilize		Do not use if package is damaged
	Use-by date		Sterilized using ethylene oxide
	Catalogue number		Batch code
	Caution		Upper limit of temperature
	Manufacturer		MR Safe

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