GYNECARE TVT™

with abdominal guides
Tension-free Support for Incontinence

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

P17170
LAB100117642v7
02/2021
Made in Switzerland
© Ethicon, Inc. 2021

Ethicon, SÀRL
Puits Godet 20
Neuchâtel
CH-2000
Switzerland
1-877-ETHICON
+1-513-337-6928

GYNECARE TVT™ with Abdominal Guides Tension-free Support for Incontinence

GYNECARE TVT Single Use Device GYNECARE TVT Reusable Introducer GYNECARE TVT Reusable Rigid Catheter Guide GYNECARE TVT Abdominal Guides and Couplers

Please read all information carefully.

Failure to properly follow instructions may result in improper functioning of the device and/or may lead to injury. Properly following instructions will not eliminate the risk of adverse reactions.

Introduction

This package insert is designed to provide instructions for use of the GYNECARE TVT with *Abdominal Guides* Tension-free Support for Incontinence System, including Single Use Device, Reusable Introducer, Reusable Rigid Catheter Guide and disposable Abdominal Guides and Couplers. It is not a comprehensive reference to surgical technique for correcting SUI (Stress Urinary Incontinence). This is not a technique manual nor a substitute for appropriate training and experience in surgical technique for correcting SUI. The device is intended to be used only by surgeons trained in the surgical treatment of SUI and specifically in implanting the GYNECARE TVT Device. Training should be in line with the training guidelines specified by the International Urogynecological Association (IUGA) or by other international, national, or local Institutions. Adverse reactions may occur notwithstanding the training and experience of the surgeon. These instructions are recommended for general use of the device. The removal of the product from the patients (if needed) should be performed by a trained urogynecologist, urologist or gynecologist experienced in the removal of pelvic floor meshes. Variations in use may occur in specific procedures due to individual technique and patient anatomy.

Description

GYNECARE TVT DEVICE

The GYNECARE TVT Tension-free Support for Incontinence device is a sterile single use device, consisting of one piece of undyed or blue (up to 0.28 weight percent phthalocyanine blue, Color Index number 74160) PROLENE™ Polypropylene Mesh (tape) approximately 1/2 inch x 18 inches (1.1 cm x 45 cm), covered by clear plastic Sheaths overlapping in the middle, and held between two stainless steel needles bonded to the mesh and Sheath with plastic collars.

Based upon the maximum implanted mesh length, for the blue mesh, the maximum amount of dye (phthalocyanine blue, Color Index number 74160) is 0.00077 g.

The GYNECARE TVT Device is available in either mechanical cut or laser cut versions for the surgeon's preference. To determine if the GYNECARE TVT Device implant is mechanical or laser cut, consult the product code on the device packaging; an (L) at the end of the number indicates the laser cut mesh.

PROLENE Polypropylene Mesh is constructed of knitted filaments of extruded polypropylene strands identical in composition to that used in PROLENE Polypropylene nonabsorbable surgical suture. The mesh is approximately 0.027 inches (0.7 mm) thick. This material, when used as a suture, has been reported to be non-reactive and to retain its strength indefinitely in clinical use. PROLENE Mesh is knitted by a process which interlinks each fiber junction.

GYNECARE TVT INTRODUCER

(Available separately – not included in GYNECARE TVT Device)

The GYNECARE TVT Introducer is provided non-sterile and is reusable. The Introducer is made of stainless steel. It consists of two parts, a handle and an inserted threaded metal shaft. The Introducer is intended to facilitate the passage of the GYNECARE TVT Device from the vagina to the abdominal skin. It is connected and fixed to the needle, via the threaded end of the shaft, prior to inserting the needle with the tape.

The GYNECARE TVT Introducer must be cleaned and sterilized prior to use. Instructions can be found later in this document.

The cleaned and sterilized GYNECARE TVT Introducer must be assembled at the point-of-use.

GYNECARE TVT RIGID CATHETER GUIDE

(Available separately – not included in GYNECARE TVT Device)

The GYNECARE TVT Rigid Catheter Guide is a non-sterile reusable instrument intended to facilitate the identification of the urethra and the bladder neck during the surgical procedure. It is inserted into a Foley catheter (recommended size 18 French) positioned in the bladder via the urethra. To facilitate

insertion, it can be lubricated with gel.

The Rigid Catheter Guide must be cleaned and sterilized prior to use. Instructions can be found in the Instructions for Cleaning Reusable Instruments Section.

GYNECARE TVT ABDOMINAL GUIDE

The GYNECARE TVT Abdominal Guide is a sterile disposable instrument intended to facilitate passage of the GYNECARE TVT Device. Two Abdominal Guides are included in each kit with the GYNECARE TVT Couplers.

GYNECARE TVT COUPLER

The GYNECARE TVT Coupler is a sterile disposable polypropylene connector used to connect the GYNECARE TVT Abdominal Guide to the GYNECARE TVT Needle. Two Couplers are included in each kit with the Abdominal Guides.

Indications / Intended Use

The GYNECARE TVT Device is intended to be used as a pubo-urethral sling for treatment of stress urinary incontinence (SUI), for female urinary incontinence resulting from urethral hypermobility and/or intrinsic sphincter deficiency. The GYNECARE TVT Introducer and Rigid Catheter Guide are available separately and GYNECARE TVT Abdominal Guides and Couplers are included in each kit.

The GYNECARE TVT Introducer is a reusable device intended to aid in the placement of the GYNECARE TVT Device retropublically.

The GYNECARE TVT Rigid Catheter Guide is a reusable device intended to facilitate the identification of the urethra and bladder neck during the placement of the GYNECARE TVT Device.

The GYNECARE TVT Abdominal Guides and Couplers are single use devices used to facilitate placement of the GYNECARE TVT Device when placed in a top-down retropublic fashion (also known as an abdominal approach).

Contraindications

As with any suspension surgery, this procedure should not be performed in pregnant patients. Additionally, because the PROLENE Polypropylene Mesh will not stretch significantly, this procedure should not be performed in patients with future growth potential, including women with plans for future pregnancy.

Patient Target Group and Factors

The GYNECARE TVT Device is intended for use for non-pregnant female patients not planning future pregnancies, who are affected by stress urinary incontinence resulting from urethral hypermobility and/or intrinsic sphincter deficiency. Surgeons should use their surgical experience and judgment to determine if PROLENE Mesh is appropriate for certain patients. Patient-specific factors may impair wound healing, which may increase the likelihood of adverse reactions.

Application / Instructions for Use

The procedure can be carried out under local anesthesia, but it can also be performed using regional or general anesthesia. The patient should be placed in the lithotomy position, taking care to avoid hip flexion greater than 60°. Insert an 18 French Foley catheter into the bladder and leave it to open drainage. The extent of dissection is minimal, i.e., a vaginal midline entry with a small para-urethral dissection to initially position the needle and two supra-pubic skin incisions. Using forceps, grasp the vaginal wall at each side of the urethra. Using a small scalpel, make a sagittal incision about 1.5 cm long starting approximately 1.0 cm from the outer urethral meatus. This incision will cover the mid urethral zone and will allow for subsequent passage of the sling (tape). With a small pair of blunt scissors, two small para-urethral dissections (approximately 0.5 cm) are made so that the tip of the needle can then be introduced into or passed through the para-urethral dissection. Then, two abdominal skin incisions of 0.5 cm –1 cm are made, one on each side of the midline, just above the symphysis, not more than 4 cm –5 cm apart. Incision placement and Needle passage near the midline and close to the back of the pubic bone are important to avoid anatomic structures in the inguinal area and lateral pelvic sidewall.

The GYNECARE TVT Rigid Catheter Guide is inserted into the channel of the Foley catheter (18 French). The handle of the Guide is fixed around the catheter, proximal to its widening. The purpose of the Guide is to move the bladder neck and urethra away from where the tip of the Needle will pass into the retropublic space by gently pushing the tip of the Foley catheter toward the posterior lateral wall of the bladder opposite to the intended Needle passage. For example, by pushing toward the patient's left side, the bladder will go from a spherical to a spheroid configuration. This moves the bladder away from the back of the pubic symphysis. Additionally, it moves the bladder neck and the urethra to the left as the Needle is passed on the patient's right side, thereby minimizing the risk of bladder perforation. The Rigid Catheter Guide, the urethra and bladder are moved contralaterally to the side of the Needle passage. During this maneuver, the bladder should be empty.

The threaded end of the Introducer is screwed into the end of one of the Needles. If the abdominal access device is utilized, the tip of the GYNECARE TVT Needle is inserted firmly into the distal (wider) end of the Coupler. Ensure there is a snug connection between the Coupler and GYNECARE TVT Needle.

Abdominal Approach

GYNECARE TVT Tension-Free Vaginal Tape can be placed in position by either an abdominal or a vaginal approach. If the abdominal approach is utilized, follow instructions above for para-urethral dissection. Holding the Abdominal Guide, gently insert it in a vertical, downward and sagittal direction, through one of the abdominal incisions until the rectus fascia is reached. Firm pressure on the Abdominal Guide will pierce the rectus fascia and bring the Abdominal Guide into the retropubic space. Upon entering the retropubic space, the operator will feel a significant decrease in resistance. Pass the Abdominal Guide through the retropubic space taking care to keep it close to the posterior surface of the pubic bone. Using a vaginal finger (of the non-dominant hand), palpate the descent of the Abdominal Guide along the posterior surface of the pubic bone. When the tip of the Abdominal Guide is palpable just inferior to the pubic bone and superior to the anterior vaginal wall, rotate the tip of the Abdominal Guide towards the midline. Confirm that you are 1 cm lateral to the Catheter Guide (mid-urethra). Using digital counter-traction against the urogenital diaphragm, guide it out through the diaphragm and out the midline vaginal incision. Once the Abdominal Guide is brought through the vaginal incision, cystoscopy should be performed to ascertain that no bladder or urethral injury has taken place.

Empty the bladder once bladder integrity has been confirmed. Screw the threaded end of the Introducer into the end of one of the Needles. Insert the tip of the Needle firmly into the distal (wider) end of the Coupler. Next pick up the GYNECARE TVT Device with attached Coupler, and holding the Abdominal Guide steady, firmly insert the tapered end of the Coupler (with attached Needle) over the Guide. Ensure there is a snug connection between the Coupler and Abdominal Guide. Be sure that this interlocked system (the Abdominal Guide, GYNECARE TVT Needle, Coupler and reusable Introducer) is oriented in the same plane. Using the non-dominant hand on the abdominal portion of the Guide and the dominant hand gripping the Introducer, **push** the system upward toward the abdomen following the curvature of the interlocked system. **The Abdominal Guide should not be used to pull the interlocked system upward toward the abdomen.** The interlocked system should be pushed up toward the abdomen until the tip of the GYNECARE TVT Needle is visible through the abdominal incision. Disarticulate the reusable Introducer and pull the remaining portion of the GYNECARE TVT Needle through the abdominal incision. Cut the tape close to the Needles. The exact same procedure is carried out on the other side.

Pull the abdominal ends of the tape upward to bring the vaginal aspect of the tape (sling) loosely, i.e., without tension, under the mid-urethra. Now, adjust the tape so that leakage is reduced allowing a few drops of urinary leakage to occur under stress. For this, use patient feedback, i.e., coughing with a full bladder (approximately 300 mL). In case of general anesthesia, verify tension of the tape exerting manual pressure on the suprapubic area with a full bladder (approximately 300 mL). Keep the vaginal incision temporarily closed by a gentle grip with small forceps. The plastic Implant Sheaths that surround the tape are then removed. **To avoid putting tension on the tape, a blunt instrument (scissors or forceps) should be placed between the urethra and the tape during removal of the plastic Implant Sheaths.** Premature removal of the Sheath may make subsequent adjustments difficult. After proper adjustment of the tape, close the vaginal incision. Cut the abdominal ends of the tape so that the ends are in the subcutis. Do not suture them. Close the skin incisions with suture or surgical skin adhesive. Empty the bladder. Following this procedure, post-operative catheterization is not typically required. The patient should be encouraged to try to empty the bladder 2–3 hours after the operation.

Vaginal Approach

If the vaginal approach is used for placement of the GYNECARE TVT Tension-Free Vaginal Tape, the Abdominal Guides and Couplers are not necessary.

Using the Introducer, the needle is passed para-urethrally, penetrating the urogenital diaphragm. Insertion and passage are controlled by using the long or index finger in the vagina under the vaginal wall on the ipsilateral side and fingertip control on the pelvic rim. The curved part of the Needle should rest in the palm of the “vaginal” hand. If you are right handed, this means that the left hand generally is the one to be used for needle guidance. With the other hand, grip the handle of the Introducer gently. Now introduce the Needle tip into the retropubic space. Once again, observe that this should be done by the palm of the vaginal hand and with the Needle tip horizontally, i.e., in the frontal plane. After passage of the urogenital diaphragm, you will feel that the resistance is significantly reduced.

Immediately aim the tip of the Needle towards the abdominal midline and lower the handle of the Introducer, thereby pressing the tip of the needle against the back of the pubic bone. Now, move the Needle tip upwards to the abdominal skin incision, keeping in close contact with the pubic bone

all the way. When the Needle tip has reached the abdominal incision, cystoscopy is performed to confirm bladder integrity. The bladder must be emptied after the first cystoscopy. Disarticulate the reusable Introducer and pull the remaining portion of the GYNECARE TVT Needle through the abdominal incision. The procedure is then repeated on the other side. The Needles are then pulled upward to bring the tape (sling) loosely, i.e., without tension, under the mid-urethra. Cut the tape close to the Needles. Now, adjust the tape so that leakage is reduced allowing a few drops of urinary leakage to occur under stress. For this, use patient feedback, i.e., coughing with a full bladder (approximately 300 mL) and keep the vaginal incision temporarily closed by a gentle grip with small forceps. The plastic Implant Sheaths that surround the tape are then removed. To avoid putting tension on the tape, a blunt instrument (scissors or forceps) should be placed between the urethra and the tape during removal of the plastic Implant Sheaths. Premature removal of the Sheath may make subsequent adjustments difficult. After proper adjustment of the tape, close the vaginal incision. The abdominal ends of the tape are then cut and left in subcutis. Do not suture them. Suture the skin incisions. Empty the bladder. Following this procedure, post-operative catheterization is not typically required. The patient should be encouraged to try to empty the bladder 2–3 hours after the operation.

Disposal

Discard unintentionally opened/ partially used/ used devices and packaging according to your facility's policies and standard operating procedures concerning biohazardous materials and waste.

Warnings and Precautions

- Do not use GYNECARE TVT Device in patients who are on anti-coagulation therapy.
- Do not use GYNECARE TVT Device in a patient who has a urinary tract infection.
- Surgeons should be familiar with surgical technique for mid-urethral sling placement and should be adequately trained in the GYNECARE TVT Device implantation procedure before employing the GYNECARE TVT Device. It is important that the Implant be located without tension under mid-urethra.
- The Abdominal Guide should not be used to pull the interlocked system upward toward the abdomen.
- Ensure there is a snug connection between the Guide and Coupler and the Coupler and GYNECARE TVT Needle.
- Acceptable surgical practice should be followed for the procedure as well as for the management of contaminated or infected wounds.
- The procedure should be performed with care to avoid large vessels, nerves, bladder and bowel. Attention to local anatomy and proper passage of needles may minimize certain risks.
- Retropubic bleeding may occur post-operatively. Observe for any symptoms or signs before releasing the patient from hospital.
- Cystoscopy should be performed to confirm bladder integrity or recognize a bladder perforation.
- The Rigid Catheter Guide should then be gently pushed into the Foley catheter so that the Catheter Guide does not extend into the holes of the Foley Catheter.
- When removing the Rigid Catheter Guide, open the handle completely so that the catheter remains properly in place.
- Do not remove the plastic Sheath until the Implant has been properly positioned.
- Ensure that the Implant is placed with minimal tension under mid-urethra.
- Do not perform this procedure if you think the surgical site may be infected or contaminated.
- Since limited clinical information is available about pregnancy following pubo-urethral sling procedure with the GYNECARE TVT Device, the patient should be counseled that future pregnancies may negate the effects of the surgical procedure and the patient may again become incontinent.
- Since limited clinical information is available with delivery following a pubo-urethral procedure with the GYNECARE TVT Device, in case of pregnancy, mode of delivery should be determined by the obstetrician in consultation with the surgeon.
- Post-operatively the patient is recommended to refrain from heavy lifting and/or exercise (e.g., cycling, jogging) for at least four to six weeks and intercourse for at least one month. The patient can return to other normal activity after one or two weeks.
- Should dysuria, bleeding or other problems occur, the patient is instructed to contact the surgeon immediately.
- All surgical instruments are subject to wear and damage under normal use. Before use, the instrument should be visually inspected. Defective instruments or instruments that appear to be corroded should not be used and should be discarded.
- As with other incontinence procedures, de novo detrusor instability may occur following a

pubo-urethral procedure utilizing the GYNECARE TVT Device. To minimize this risk, make sure to place the Implant tension-free in the mid-urethral position, as described above.

- Do not contact the PROLENE Mesh with any staples, clips or clamps as mechanical damage to the mesh may occur.
- Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and surgeons.
- Prophylactic antibiotics can be administered according to the surgeon's usual practice.
- In the event of a damaged product, the product should be discarded and a new one retrieved to complete the procedure. In the event of a product malfunction during use, continue or discontinue use at surgeon's discretion.

Adverse Reactions / Undesirable Side Effects

- Punctures or lacerations of vessels, nerves, structures, or organs, including the bladder, urethra, or bowel, may occur and may require surgical repair.
- Local irritation at the wound site may occur.
- As with any Implant, a foreign body response will occur, the extent of which may differ. This response could result in extrusion, erosion, exposure, fistula formation and/or inflammation, the severity of which is unpredictable, or other adverse reactions, which may be ongoing.
- Fistula formation, acute and chronic inflammation and ongoing risk of mesh extrusion, exposure, or erosion into the vagina or other structures or organs (such as bladder, urethra or rectum) which may be difficult to treat and result in consequent pain. Mesh extrusion, exposure, or erosion into the vagina may also cause offensive vaginal discharge.
- Infection following transvaginal implantation. As with all surgical procedures and the implantation of foreign bodies, there is a risk of infection and PROLENE Mesh may potentiate an existing infection.
- Pain - which may be severe and chronic.
- Temporary or chronic voiding dysfunction (or difficulty voiding) or urinary retention/obstruction independent from that caused by overcorrection or urethral hypermobility, i.e. too much tension applied to the tape, or from misplacement of the sling or placing the sling too tightly.
- Pain with intercourse, (dyspareunia) and loss of sexual function, (apareunia), which may be ongoing and may not resolve in some patients.
- Excessive contraction or shrinkage of the tissue surrounding the mesh, and vaginal scarring from causes which include, but are not limited to, chronic inflammation and mesh exposure.
- Neuromuscular problems, including acute and/or chronic pain in the groin, thigh, leg, pelvic and/or abdominal area, and leg weakness, may occur.
- Recurrence of incontinence.
- Bleeding including hemorrhage, or hematoma.
- One or more revision surgeries may be necessary to treat these adverse reactions. Revision surgeries may not resolve complications and are associated with a risk of adverse reactions.
- PROLENE Mesh is a permanent implant that integrates into the tissue. In cases in which the PROLENE Mesh needs to be removed in part or whole, significant dissection may be required.

Removal of the Implant in whole or in part will not necessarily alleviate the patient's symptoms. Removal of part of the Implant can be difficult. Surgery to remove the whole or part of an Implant can result in further scarring and tissue damage which, in turn, may have adverse outcomes including severe chronic pain which may not be able to be satisfactorily treated. Surgery to remove the whole or part of the Implant may also result in recurrence of SUI. Removal of the eroded mesh will not necessarily prevent further erosions or other adverse events.

Other Adverse Reaction / Undesirable Side Effects

- Seroma
- Urge incontinence, including de novo urge incontinence
- Urinary frequency
- Urinary retention
- Adhesion formation
- Atypical vaginal discharge
- Exposed mesh may cause pain or discomfort to the patient's partner during intercourse.
- Death

The surgeon should convey adverse reactions, undesirable side effects and risks associated with the product and the procedure to the patient prior to undertaking SUI and advise the patient to contact a surgeon in case of any deviation from the normal post-operative course. Any serious incident that

has occurred in relation to the device should be reported to the manufacturer and the country-competent authority.

Magnetic Resonance Imaging (MRI)

GYNECARE TVT Device is Magnetic Resonance (MR) Safe.

Performance(s) / Action(s)

Animal studies show that implantation of PROLENE Mesh elicits an inflammatory reaction in tissues, which is transient and is followed by the deposition of a thin fibrous layer of tissue, which can grow through the interstices of the mesh, thus incorporating the mesh into adjacent tissue. The material is not absorbed, nor is it subject to degradation or weakening by the action of tissue enzymes. The device is non-absorbable therefore it remains for the lifetime of the patient unless removed.

Sterility

The GYNECARE TVT Device is sterile if kept in original, unopened packaging. DO NOT RESTERILIZE. DO NOT REUSE. Do not use if package is opened or damaged.

The reusable GYNECARE TVT Introducer and GYNECARE TVT Rigid Catheter Guide are supplied separately and are non-sterile. They are to be cleaned and sterilized prior to each use as described below.

Instructions for Cleaning Reusable Instruments (available separately)

(GYNECARE TVT Introducer, GYNECARE TVT Rigid Catheter Guide)

To ensure the reliability and functionality of GYNECARE TVT Introducer and GYNECARE TVT Rigid Catheter Guide, clean and sterilize the instruments before initial use and after each procedure. The following are suggested pre-cleaning and automated methods for cleaning the instruments. Each instrument must be treated individually and matt surfaces must be protected against scratches. The design of this instrument permits effective cleaning and reliable sterilization. The user must ensure that cleaning and sterilization are performed in accordance with guidelines, standards and National Health Authority requirements.

Prior to cleaning, the GYNECARE TVT Introducer should be separated into its component parts (handle and threaded shaft). The Introducer is reassembled after cleaning and sterilization prior to use.

Point-of-Use Care of the Instruments

Disassemble the instrument and wipe clean of visible soil. Reassemble the instrument and properly contain it for transport to a dedicated cleaning work area, ensuring the instrument is kept moist to avoid drying of soil.

Pre-Cleaning of the Instruments

1. The GYNECARE TVT Introducer should be separated into its component parts (handle and threaded shaft). The area of the Rigid Catheter Guide that clips onto the long shaft should be opened so that all surfaces are easily accessible for cleaning.
2. Rinse the instrument thoroughly with tap water.
3. Prepare a neutral pH enzymatic solution suitable for stainless steel instruments per the manufacturer's recommendations.
4. Soak the instrument components in the enzymatic cleaner for a minimum of 10 minutes.
5. Wash in the enzymatic solution at a temperature recommended by the manufacturer. Remove any contamination from body fluids or tissues using a soft-bristled brush (Spectrum M16 or equivalent) and scrubbing the exposed surfaces of the GYNECARE TVT Introducer handle and threaded shaft for a minimum of 2 minutes. The lumen of the handle should be cleaned with a round soft-bristled brush by keeping the instrument submerged in the cleaning solution and scrubbing for at least 2 minutes. The Rigid Catheter Guide should be cleaned using a soft-bristled brush (Spectrum M16 or equivalent) and scrubbing the surfaces for a minimum of 2 minutes.
6. Remove the instrument from the enzymatic solution and rinse thoroughly under running tap water for a minimum of 1 minute.
7. Place the instrument in an ultrasonic cleaner with fresh enzymatic detergent for a minimum of 10 minutes. The frequency should be 25 kHz to 50 kHz (target 35 kHz).
8. Remove the instrument from the ultrasonic cleaner and rinse thoroughly under running tap water for a minimum of 1 minute.
9. Visually inspect the instrument for cleanliness. If necessary, repeat steps 1-8. Following pre-cleaning, the instrument should be cleaned in a washer-disinfector that conforms to ISO 15883-1 and 2.

Automated Cleaning Method:

After the pre-cleaning steps have been carried out, the instruments can be cleaned in an automated washing cycle. Automatic washing cycles are suitable for stainless steel instruments. One recommended cycle is described below:

PHASE	MINIMUM TIME (MINUTES)	TEMPERATURE AND WATER QUALITY	DETERGENT TYPE AND CONCENTRATION
Pre-wash	02:00	Cold tap water	N/A
Enzyme wash	15:00	Hot tap water or detergent manufacturer's recommendation Minimum 176°F (80°C)	Enzymatic Detergent (Enzol 2 oz./gallon or an equivalent product)
Rinse 1	02:00	Hot tap water	N/A
Rinse 2	15:00	Hot tap water	N/A
Final rinse	05:00	Hot tap water	N/A
Rinse with demineralized or purified water	02:00	Demineralized or purified water 176°F (80°C)	N/A
Thermal disinfection (optional)	01:00	Demineralized or purified water 194°F (90°C)	N/A
Drying	10:00-15:00	199.4°F (93°C)	N/A

Thermal Disinfection (Optional)

Note that disinfection is mandatory in some areas of the world, e.g., EU. The user must ensure that cleaning, disinfection and sterilization are conducted in accordance with the appropriate guidelines, standards, or National Health Authority requirements.

Thermal disinfection can be included as part of the automatic cleaning cycle. The instrument can be thermally disinfected at 194°F (90°C) or greater for a minimum of 1 minute to render the instrument safe for handling. The GYNECARE TVT Introducer should be disassembled for thermal disinfection.

Note: Do not exceed 212°F (100°C) for longer than 10 minutes.

**Sterilization Recommendations for Reusable Instruments (available separately)
(GYNECARE TVT Introducer, GYNECARE TVT Rigid Catheter Guide)**

The GYNECARE TVT Introducer and GYNECARE TVT Rigid Catheter Guide are supplied non-sterile. To sterilize, steam autoclave prior to each use. Steam (moist heat) sterilization is recommended in an approved, pre-vacuum cycle listed below. The steam sterilizer must be validated to the requirements of any local standards and guidance such as EN285 or AAMI/ANSI ST8. Ensure that the steam sterilizer cycle is chosen that is designed to remove air from porous loads and/or lumened instruments in accordance to manufacturer's instructions.

Prior to steam sterilization ensure the instrument is dry and package the instrument in two layers of sterilization wrap or other suitable sterile barrier package. The sterile barrier packaging should meet the requirements in ISO 11607-1.

WARNING: The GYNECARE TVT Introducer must be sterilized disassembled.

WARNING: The section of the Rigid Catheter Guide that clips onto the long shaft must be open during cleaning and sterilization. If it is not open during cleaning and sterilization, the cleaning process may be less effective, and sterility of the instrument may be compromised.

The types of sterilization equipment and available cycles may vary by geographic region. The following sterilization parameters have been validated to deliver a sterility assurance level of 10⁻⁶ for this instrument:

Market Area	Cycle Type	Temperature Set Point	Exposure Time (minutes)	Minimum Drying Time (minutes)
USA and its territories	Pre-vacuum	270°F/132°C	4	35
EU and others	Pre-vacuum	273.2°F/134°C	3	35

Note:

The above table includes the minimum temperature and cycle time validated to assure sterility.

It is the responsibility of the end user to assure sterility of the product when using sterilization process recommended, since bioburden and sterilization equipment will vary.

Instrument Maintenance

- GYNECARE TVT Introducer
Before each use, inspect the threaded parts of the inner shaft.
- GYNECARE TVT Rigid Catheter Guide
Before each use, inspect the instrument. Check to ensure that the long end, which traverses the catheter channel, has no sharp edges or burrs.

How Supplied

The GYNECARE TVT Device and Abdominal Guides and Couplers are provided sterile (ethylene oxide) for single use. The reusable GYNECARE TVT Introducer and GYNECARE TVT Rigid Catheter Guide are supplied separately and are non-sterile. These accessories are to be cleaned and sterilized prior to each use as described above.

Storage

No special storage conditions required. Do not use after expiration date.

Traceability

A traceability label that identifies the type, size, expiry date, and batch code of the device is included in every package. This label should be affixed or electronically added to the patient’s permanent medical record to clearly identify the device that was implanted.

Symbols Used on Labeling

	Caution		Catalogue number
	Do not re-use		Sterilized using ethylene oxide
	Do not resterilize		Manufacturer
	Do not use if package is damaged		MR Safe
	Use-by date		Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	Batch code		

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