

COMPETITIVE COMPARISON¹

	Caldera Medical Desara[®] One	Coloplast Altis[®]	Boston Scientific Solyx[™]
One handed deployment	✓		
Anchor Anti-rotation mechanism	✓		
One handed anchor release mechanism	✓		
Large Pore Size ($\approx > 1000 \mu\text{m}$) and high porosity	✓		✓
Blue Mesh	✓		✓
Fully de-tanged edges	✓		

PRODUCT REFERENCE

ORDER NUMBER	DESCRIPTION	QUANTITY
CAL-DS2113	Desara [®] One Single Incision Sling System for Female Stress Urinary Incontinence	Box of 3



For product evaluation, call your Caldera Sales Representative today at (866) 4-CALDERA




For more information visit: www.calderamedical.com



Contact us at:
Phone 818.879.6555
Fax 818.879.6556




The **Single Incision Sling**
you asked for



Bidirectional adjustability

Disengagement wings on introducer provide anchor continuity for easy tensioning prior deployment.¹



Anchor anti-rotation

Square shape insertion pin, minimizes potential for mesh twisting during delivery.¹



Midline mark



One handed anchor deployment system



Unique Introducer

The ergonomic and lightweight design allows for single handed deployment.^{1,2}


Desara[®] One is the choice for single incision slings*, providing simplified tensioning and streamlined implantation.¹



High pull-out force¹



Desara[®] Blue large pore size mesh



Desara[®] Blue mesh test:¹
Before and after of 5 lbs of force applied.



Anchor fully integrated with mesh

Helps prevent mesh ripping, twisting and roping, maximizing the strength of the connection.¹



Streamlined anchor-mesh transition

Designed to minimize dissection, lessen tissue disruption and to lay flat under the urethra.¹



LITERATURE

- 1. Caldera Medical data on file.
- 2. Desara[®] One IFU.

*Single incision slings may have a lower cure rate compared with standard length mid-urethral slings.