


# Benesta®

HYSTEROSCOPIC TISSUE REMOVAL SYSTEM

Battery powered system puts control at your fingertips




Intuitive **push to cut** trigger allows for excellent control<sup>1</sup>

Single use **battery powered** system

1mm close-cutting edge for improved access to **hard-to-reach** fundal and cornual tissue

15mm



Long working length provides better and comfortable access for high BMI patients<sup>1</sup>


228mm

Optimized for use with a **hanging bag**

Small 5.8mm **High Definition**, 0° rod lens scope with beveled distal tip for comfortable insertion

 1. Caldera Medical Data on File

## SYSTEM SET-UP



## ORDER INFORMATION

ORDER NUMBER	PRODUCT DESCRIPTION	QUANTITY
Benesta® Hysteroscope Set and Sterilization Tray		
CAL-TR111	Benesta® Hysteroscope Set (includes outflow channel, light post adapters and cleaning brushes)	1 Each
CAL-TR131	Benesta® Sterilization Tray	Single Unit
Benesta® Consumables		
CAL-TR1513	Benesta® Tissue Removal Device (pack of 3)	3 pack
CAL-TR1415	Single Use Fluid Management pack (1 Y-tubing + 2 Single Use seals)	5 pack
CAL-14210	Single use seals only, pack of 10	10 pack



For product evaluation, call your Caldera Sales Representative today at (866) 4-CALDERA



For more information visit:  
[www.calderamedical.com](http://www.calderamedical.com)



Contact us at:  
Phone 818.879.6555  
Fax 818.879.6556